

HuC/HuD Rabbit mAb


ABZ2034

Product Overview

Name	HuC/HuD Rabbit mAb
Accession(Primary)	Q60900
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 40 kDa; Observed MW: 40 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:500-1:2000; ICC/IF-1:200-1:500; FC-1:200-1:500; IP-1:20-1:50

Target information (Q60900)

Gene Name	Elavl3/Elavl4
About Gene Symbol	



Western blot analysis of HuD/HuC in rat Brain
lysates using HuC/HuD antibody.