

KLF4 Rabbit mAb


ABZ2035

Product Overview


Name	KLF4 Rabbit mAb
Accession(Primary)	Q60793
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Mouse, Rat
Calculated MW	Calculated MW: 52 kDa; Observed MW: 52 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:1000; ICC/IF-1:500-1:1000; IP-1:20-1:50

Target information (Q60793)

Gene Name	Klf4
About Gene Symbol	



Immunohistochemistry analysis of paraffin-embedded mouse colon using KLF4 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of KLF4 in HeLa, CHO-K1 lysates using KLF4 antibody.