

NDUFAF5 Rabbit mAb


ABZ2036

Product Overview


Name	NDUFAF5 Rabbit mAb
Accession(Primary)	Q5TEU4
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 39 kDa; Observed MW: 39 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:1000; ICC/IF-1:200-1:500; FC-1:100-1:200; IP-1:50-1:100

Target information (Q5TEU4)

Gene Name	NDUFAF5
About Gene Symbol	



Western blot analysis of NDUF5 in 3T3 lysates using NDUF5 antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using NDUF5 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.