

Artemin Rabbit mAb

ABZ2037


Product Overview

Name	Artemin Rabbit mAb
Accession(Primary)	Q5T4W7
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 23 kDa; Observed MW: 23 kDa
Recommended Dilution	WB-1:1000-1:5000; IP-1:10-1:100; IHC-1:50-1:100


Target information (Q5T4W7)

Gene Name ARTN

About Gene Symbol



Western blot analysis of Artemin in K562, rat Brain, C6, 3T3, HeLa lysates using Artemin antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using Artemin antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.