

PDZK1 Rabbit mAb

ABZ2038


Product Overview

Name	PDZK1 Rabbit mAb
Accession(Primary)	Q5T2W1
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 57 kDa; Observed MW: 70 kDa
Recommended Dilution	WB-1:1000-1:5000; IP-1:20-1:50

Target information (Q5T2W1)

Gene Name PDZK1

About Gene Symbol



Western blot analysis of PDZK1 in T4D7, Caco2 lysates using PDZK1 antibody.

