


DPS1 Rabbit mAb

ABZ2039

Product Overview

Name	DPS1 Rabbit mAb
Accession(Primary)	Q5T2R2
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 46 kDa; Observed MW: 46 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20-1:50

Target information (Q5T2R2)

Gene Name PDSS1**About Gene Symbol**

Western blot analysis of DPS1 in K562 lysates
using DPS1 antibody.

