

Transmembrane Protein 85 Rabbit mAb


ABZ2041

Product Overview


Name	Transmembrane Protein 85 Rabbit mAb
Accession(Primary)	Q5J8M3
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 20 kDa; Observed MW: 20 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; IP-1:50-1:100

Target information (Q5J8M3)

Gene Name	EMC4
About Gene Symbol	



Western blot analysis of TMEM85 in 3T3 lysates using Transmembrane Protein 85 antibody.



Western blot analysis of TMEM85 in Jurkat, rat Brain, C6, CHO-K1, HeLa lysates using TMEM85 antibody.