

Gli3 Rabbit mAb


ABZ2042

Product Overview

Name	Gli3 Rabbit mAb
Accession(Primary)	Q5IS56
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 170 kDa; Observed MW: 200 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:50-1:100; ICC/IF-1:50-1:100

Target information (Q5IS56)

Gene Name	GLI3
About Gene Symbol	



Western blot analysis of Gli3 in 293 lysates using
GLI3 antibody.