


MTX3 Rabbit mAb


ABZ2043

Product Overview

Name	MTX3 Rabbit mAb
Accession(Primary)	Q5HYI7
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 35 kDa; Observed MW: 35 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (Q5HYI7)

Gene Name MTX3**About Gene Symbol**



Western blot analysis of MTX3 in mouse brain
lysates using MTX3 antibody.