

## Meckelin Rabbit mAb


ABZ2044

### Product Overview


|                             |                                                                            |
|-----------------------------|----------------------------------------------------------------------------|
| <b>Name</b>                 | Meckelin Rabbit mAb                                                        |
| <b>Accession(Primary)</b>   | Q5HYA8                                                                     |
| <b>Host Species</b>         | Rabbit                                                                     |
| <b>Clonality</b>            | Monoclonal                                                                 |
| <b>Isotype</b>              | IgG                                                                        |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                                                          |
| <b>Calculated MW</b>        | Calculated MW: 112 kDa; Observed MW: 112 kDa                               |
| <b>Recommended Dilution</b> | WB-1:500-1:1000; IHC-P-1:100-1:200;<br>ICC/IF-1:500-1:1000; FC-1:200-1:500 |

### Target information (Q5HYA8)

|                          |        |
|--------------------------|--------|
| <b>Gene Name</b>         | TMEM67 |
| <b>About Gene Symbol</b> |        |



Western blot analysis of Meckelin in Hela, A549, HL-60, U251, U87-MG lysates using Meckelin antibody.



Immunocytochemistry analysis of Meckelin (green) in MCF-7 using Meckelin antibody, and DAPI (blue)