

14-3-3 epsilon Rabbit mAb


ABZ2046

Product Overview

Name	14-3-3 epsilon Rabbit mAb
Accession(Primary)	P62258
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW:29 kDa;Observed MW: 29 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:100-1:200; FC-1:100-1:1000

Target information (P62258)

Gene Name	YWHAE
About Gene Symbol	



Western blot analysis of 14-3-3 epsilon in 293T, SW480 lysates using 14-3-3 epsilon antibody.