

BRG1 Rabbit mAb


ABZ2047

Product Overview


Name	BRG1 Rabbit mAb
Accession(Primary)	Q3TKT4
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 181 kDa; Observed MW: 220 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:200-1:500; FC-1:100-1:200; IP-1:10-1:100

Target information (Q3TKT4)

Gene Name	Smarca4
About Gene Symbol	



Western blot analysis of BRG1 in 3T3, Raw264.7, C6 lysates using BRG1 antibody.



Western blot analysis of BRG1 in 293T, HeLa, HepG2, A549 lysates using BRG1 antibody.