

HKDC1 Rabbit mAb

ABZ2048


Product Overview

Name	HKDC1 Rabbit mAb
Accession(Primary)	Q2TB90
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 103 kDa; Observed MW: 103 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20-1:50

Target information (Q2TB90)

Gene Name HKDC1

About Gene Symbol



Western blot analysis of HKDC1 in HeLa, CHO-K1, Ramos lysates using HKDC1 antibody.

