


GEN1 Rabbit mAb

ABZ2049

Product Overview

Name	GEN1 Rabbit mAb
Accession(Primary)	Q17RS7
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Rat
Calculated MW	Calculated MW: 103 kDa; Observed MW: 103 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (Q17RS7)

Gene Name GEN1**About Gene Symbol**

Western blot analysis of GEN1 in rat Brain, C6, 3T3, HeLa lysates using GEN1 antibody.

