

Mitofilin Rabbit mAb


ABZ2050

Product Overview

Name	Mitofilin Rabbit mAb
Accession(Primary)	Q16891
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 84 kDa; Observed MW: 84 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:100-1:500; FC-1:10-1:100

Target information (Q16891)

Gene Name	IMMT
About Gene Symbol	IMMT (Inner Membrane Mitochondrial Protein) is a membrane receptor involved in receptor-mediated signaling. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms an



Western blot analysis of Mitofilin in K562, Molt4, Raji lysates using Mitofilin antibody.