

TXNRD1 Rabbit mAb


ABZ2051

Product Overview


Name	TXNRD1 Rabbit mAb
Accession(Primary)	Q16881
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 71 kDa; Observed MW: 55 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:100-1:200; FC-1:100-1:500; IP-1:10-1:100

Target information (Q16881)


Gene Name	TXNRD1
About Gene Symbol	




Western blot analysis of TXNRD1 in 3T3 lysates using Thioredoxin Reductase 1 antibody.



Western blot analysis of TXNRD1 in K562, HeLa lysates using TXNRD1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using TXNRD1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of TXNRD1(green) in Hela using TXNRD1 antibody, and DAPI(blue).

