

PFKFB3 Rabbit mAb


ABZ2052

Product Overview


Name	PFKFB3 Rabbit mAb
Accession(Primary)	Q16875
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 60 kDa; Observed MW: 60 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC:1:20-1:50; ICC/IF-1:50-1:100; FC-1:20; IP-1:20-1:50

Target information (Q16875)

Gene Name	PFKFB3
About Gene Symbol	PFKFB3 (6-Phosphofructo-2-Kinase/Fructose-2,6-Biphosphatase 3) is a protein phosphatase involved in adaptive immunity. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic ta



Western blot analysis of PFKFB3 in Raw264.7, C6 lysates using PFKFB3 antibody.



Western blot analysis of PFKFB3 in K562 lysates using PFKFB3 antibody.