

DUSP6 Rabbit mAb


ABZ2053

Product Overview

Name	DUSP6 Rabbit mAb
Accession(Primary)	Q16828
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 42 kDa; Observed MW: 42 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; FC-1:50-1:200; IP-1:20-1:50


Target information (Q16828)

Gene Name	DUSP6
About Gene Symbol	




Copyright© 2024 Zen-bioscience Limited.


Western blot analysis of DUSP6 in NIH/3T3 lysates using DUSP6 antibody.



Western blot analysis of DUSP6 in NIH/3T3 lysates using DUSP6 antibody.



Immunohistochemistry analysis of paraffin-embedded Human thyroid using DUSP6 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of DUSP6 in Mouse brain
□ Rat brain lysates using DUSP6 antibody.