

## PCK2 Rabbit mAb


ABZ2054

### Product Overview



|                             |                                                                          |
|-----------------------------|--------------------------------------------------------------------------|
| <b>Name</b>                 | PCK2 Rabbit mAb                                                          |
| <b>Accession(Primary)</b>   | Q16822                                                                   |
| <b>Host Species</b>         | Rabbit                                                                   |
| <b>Clonality</b>            | Monoclonal                                                               |
| <b>Isotype</b>              | IgG                                                                      |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                                                        |
| <b>Calculated MW</b>        | Calculated MW: 71 kDa; Observed MW: 71 kDa                               |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000; IHC-P-1:200-1:500;<br>ICC/IF-1:200-1:500; FC-1:20-1:50 |

### Target information (Q16822)


|                          |      |
|--------------------------|------|
| <b>Gene Name</b>         | PCK2 |
| <b>About Gene Symbol</b> |      |



Immunocytochemistry analysis of PCK2 (green) in HT-1080 using PCK2 antibody, and DAPI (blue).



Western blot analysis of PCK2 in C6, 3T3, HeLa lysates using PCK2 antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer tissue using PCK2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.