

GSTA3 Rabbit mAb

ABZ2055


Product Overview

Name	GSTA3 Rabbit mAb
Accession(Primary)	Q16772
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 25 kDa; Observed MW: 25 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (Q16772)

Gene Name GSTA3

About Gene Symbol



Western blot analysis of GSTA3 in C6, 3T3, HeLa lysates using GSTA3 antibody.

