

CYP1B1 Rabbit mAb

ABZ2057


Product Overview

Name	CYP1B1 Rabbit mAb
Accession(Primary)	Q16678
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 61 kDa; Observed MW: 61 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500


Target information (Q16678)

Gene Name CYP1B1

About Gene Symbol



Western blot analysis of CYP1B1 in rat Brain lysates using Cytochrome P450 1B1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using CYP1B1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.