

NRF1 Rabbit mAb


ABZ2059

Product Overview


Name	NRF1 Rabbit mAb
Accession(Primary)	Q16656
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 54 kDa; Observed MW: 68 kDa
Recommended Dilution	WB-1:1000-1:10000; IHC-1:50-1:2000; ICC/IF-1:50-1:100; FC-1:10-1:100; IP-1:10-1:100

Target information (Q16656)


Gene Name	NRF1
About Gene Symbol	



Western blot analysis of NRF1 in C6, 3T3, HeLa lysates using Nrf1 antibody.



Immunocytochemistry analysis of NRF1(green) in HeLa using NRF1 antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using NRF1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.