

MOG Rabbit mAb


ABZ2060

Product Overview


Name	MOG Rabbit mAb
Accession(Primary)	Q16653
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 28 kDa; Observed MW: 28 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:500-1:5000; ICC/IF-1:50-1:100

Target information (Q16653)

Gene Name	MOG
About Gene Symbol	



Immunocytochemistry analysis of MOG (green) in SH-SY5Y using MOG antibody, and DAPI (blue).



Western blot analysis of Myelin oligodendrocyte glycoprotein in rat Brain lysates using Myelin oligodendrocyte glycoprotein antibody.