


MAPKAPK3 Rabbit mAb

ABZ2061

Product Overview

Name	MAPKAPK3 Rabbit mAb
Accession(Primary)	Q16644
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 43 kDa; Observed MW: 43 kDa
Recommended Dilution	WB-1:1000-1:5000; FC-1:50-1:100; IP-1:20-1:50

Target information (Q16644)

Gene Name MAPKAPK3**About Gene Symbol**

Western blot analysis of MK3 in A549, HL-60
lysates using MAPKAP Kinase 3 antibody.

