

SMN/Gemin 1 Rabbit mAb


ABZ2062

Product Overview


Name	SMN/Gemin 1 Rabbit mAb
Accession(Primary)	Q16637
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 32 kDa; Observed MW: 35 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500

Target information (Q16637)

Gene Name	SMN1
About Gene Symbol	SMN1 (Survival Of Motor Neuron 1, Telomeric) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of SMN/Gemin 1 in HeLa lysates using SMN antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using Gemin 1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.