

Trk B Rabbit mAb


ABZ2064

Product Overview

Name	Trk B Rabbit mAb
Accession(Primary)	Q16620
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 92 kDa; Observed MW: 120-140 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500; ICC/IF-1:100-1:200

Target information (Q16620)

Gene Name	NTRK2
About Gene Symbol	



Western blot analysis of TrkB (pan) in rat Brain
lysates using Trk B antibody.