

# Extracellular Matrix Protein 1 Rabbit mAb


ABZ2066

## Product Overview



|                             |                                                                                         |
|-----------------------------|-----------------------------------------------------------------------------------------|
| <b>Name</b>                 | Extracellular Matrix Protein 1 Rabbit mAb                                               |
| <b>Accession(Primary)</b>   | Q16610                                                                                  |
| <b>Host Species</b>         | Rabbit                                                                                  |
| <b>Clonality</b>            | Monoclonal                                                                              |
| <b>Isotype</b>              | IgG                                                                                     |
| <b>Species Reactivity</b>   | Human, Rat                                                                              |
| <b>Calculated MW</b>        | Calculated MW: 61 kDa; Observed MW: 61 kDa                                              |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000; IHC-P-1:100-1:200;<br>ICC/IF-1:100-1:200; FC-1:20-1:50; IP-1:10-1:100 |

## Target information (Q16610)


|                          |      |
|--------------------------|------|
| <b>Gene Name</b>         | ECM1 |
| <b>About Gene Symbol</b> |      |



Immunocytochemistry analysis of Extracellular Matrix Protein 1 (green) in K562 using Extracellular Matrix Protein 1 antibody, and DAPI (blue).



Western blot analysis of Extracellular matrix protein 1 in HeLa lysates using Extracellular matrix protein 1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using Extracellular matrix protein 1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.