


MLK3 Rabbit mAb

ABZ2067

Product Overview

Name	MLK3 Rabbit mAb
Accession(Primary)	Q16584
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 93 kDa; Observed MW: 93 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:100

Target information (Q16584)

Gene Name MAP3K11**About Gene Symbol**

Western blot analysis of MLK3 in HepG2, A431 lysates using MLK3 antibody.

