

Synaptophysin Like Protein 1 Rabbit mAb


ABZ2068

Product Overview


Name	Synaptophysin Like Protein 1 Rabbit mAb
Accession(Primary)	Q16563
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 29 kDa; Observed MW: 29 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; FC-1:50-1:100

Target information (Q16563)

Gene Name	SYPL1
About Gene Symbol	SYPL1 (Synaptophysin Like 1) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of SYPL1 in K562 lysates using Synaptophysin Like Protein 1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using SYPL1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.