

PPP2R5E Rabbit mAb


ABZ2073

Product Overview


Name	PPP2R5E Rabbit mAb
Accession(Primary)	Q16537
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 55 kDa; Observed MW: 55 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500; ICC/IF-1:200-1:500; FC-1:200-1:1000; IP-1:50-1:100

Target information (Q16537)

Gene Name	PPP2R5E
About Gene Symbol	PPP2R5E (Protein Phosphatase 2 Regulatory Subunit B'Epsilon) is a protein phosphatase involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellul



Western blot analysis of PPP2R5E in C6, 3T3, mouse brain, rat brain lysates using PPP2R5E antibody.



Western blot analysis of PPP2R5E in Human placenta, Hela, 293T lysates using PPP2R5E antibody.