


HLF Rabbit mAb

ABZ2074

Product Overview

Name	HLF Rabbit mAb
Accession(Primary)	Q16534
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 33 kDa; Observed MW: 33 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (Q16534)

Gene Name HLF**About Gene Symbol**

Western blot analysis of HLF in HeLa, A549, HL-60, U87-MG lysates using HLF antibody.

