

PKN1 Rabbit mAb


ABZ2076

Product Overview

Name	PKN1 Rabbit mAb
Accession(Primary)	Q16512
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 104 kDa; Observed MW: 120 kDa
Recommended Dilution	WB-1:1000-1:5000; FC-1:10-1:100

Target information (Q16512)

Gene Name	PKN1
About Gene Symbol	PKN1 (Protein Kinase N1) is a serine/threonine protein kinase involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of PKN in C6 lysates using PKN1 antibody.