

ADRM1/ARM-1 Rabbit mAb


ABZ2077

Product Overview


Name	ADRM1/ARM-1 Rabbit mAb
Accession(Primary)	Q16186
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 42 kDa; Observed MW: 42 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; ICC/IF-1:500-1:1000; FC-1:50-1:100; IP-1:10-1:100

Target information (Q16186)


Gene Name	ADRM1
About Gene Symbol	



Western blot analysis of ADRM1 in K562, HeLa lysates using ADRM1 antibody.



Immunocytochemistry analysis of ADRM1(green) in HeLa using ADRM1 antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using ADRM1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.