

SEPT7 Rabbit mAb


ABZ2078

Product Overview


Name	SEPT7 Rabbit mAb
Accession(Primary)	Q16181
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 51 kDa; Observed MW: 51 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:50-1:100

Target information (Q16181)

Gene Name	SEPTIN7
About Gene Symbol	SEPTIN7 (Septin 7) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of SEPT7 in 3T3, rat brain lysates using Septin 7 antibody.



Western blot analysis of SEPT7 in HeLa, MCF-7, mouse brain lysates using SEPT7 antibody.