

ZYX Rabbit mAb


ABZ2079

Product Overview


Name	ZYX Rabbit mAb
Accession(Primary)	Q15942
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 61 kDa; Observed MW: 78 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:1000; ICC/IF-1:200-1:500; FC-1:50-1:100; IP-1:20-1:50

Target information (Q15942)


Gene Name	ZYX
About Gene Symbol	




Immunocytochemistry analysis of ZYX (green) in Hela using ZYX antibody, and DAPI (blue).



Western blot analysis of ZYX in PC-3, HUVEC lysates using ZYX antibody.



Western blot analysis of ZYX in C6, rat Brain lysates using ZYX antibody



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using Zyxin antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.