


Adiponectin Rabbit mAb

ABZ2081

Product Overview

Name	Adiponectin Rabbit mAb
Accession(Primary)	Q15848
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 26 kDa; Observed MW: 30 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (Q15848)

Gene Name ADIPOQ**About Gene Symbol**

Western blot analysis of Adiponectin in K562, HeLa lysates using Adiponectin antibody.

