

NEDD8 Rabbit mAb


ABZ2082

Product Overview

Name	NEDD8 Rabbit mAb
Accession(Primary)	Q15843
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 9 kDa; Observed MW: 9 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:500; ICC/IF-1:100-1:200; FC-1:200-1:500; IP-1:20-1:50

Target information (Q15843)


Gene Name	NEDD8
About Gene Symbol	



Western blot analysis of NEDD8 in Jurkat lysates using NEDD8 antibody.



Immunocytochemistry analysis of NEDD8(green) in HeLa using NEDD8 antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using NEDD8 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.