

CD226 Rabbit mAb

ABZ2085


Product Overview

Name	CD226 Rabbit mAb
Accession(Primary)	Q15762
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 39 kDa; Observed MW: 60-80 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-:1:200-1:2000

Target information (Q15762)

Gene Name CD226

About Gene Symbol



Western blot analysis of CD226 in Jurkat lysates using CD226 antibody.

