

PSD93 Rabbit mAb

ABZ2088


Product Overview

Name	PSD93 Rabbit mAb
Accession(Primary)	Q15700
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 98 kDa; Observed MW: 110 kDa
Recommended Dilution	WB-1:1000-1:5000; FC-1:100-1:500

Target information (Q15700)

Gene Name DLG2

About Gene Symbol



Western blot analysis of PSD93 in HeLa, C6, Ramos lysates using PSD93 antibody.

