

## IKB beta Rabbit mAb


ABZ2089

### Product Overview


|                             |                                                                                       |
|-----------------------------|---------------------------------------------------------------------------------------|
| <b>Name</b>                 | IKB beta Rabbit mAb                                                                   |
| <b>Accession(Primary)</b>   | Q15653                                                                                |
| <b>Host Species</b>         | Rabbit                                                                                |
| <b>Clonality</b>            | Monoclonal                                                                            |
| <b>Isotype</b>              | IgG                                                                                   |
| <b>Species Reactivity</b>   | Human                                                                                 |
| <b>Calculated MW</b>        | Calculated MW: 38 kDa; Observed MW: 48 kDa                                            |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000; IHC-P-1:200-1:500;<br>ICC/IF-1:20-1:50; FC-1:20-1:50; IP-1:10-1:100 |

### Target information (Q15653)

|                          |        |
|--------------------------|--------|
| <b>Gene Name</b>         | NFKB1B |
| <b>About Gene Symbol</b> |        |



Western blot analysis of IKB beta in K562, 3T3, HeLa lysates using IKB beta antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using IKB beta antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.