

Sonic Hedgehog Rabbit mAb


ABZ2091

Product Overview

Name	Sonic Hedgehog Rabbit mAb
Accession(Primary)	Q15465
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 50 kDa; Observed MW: 50,27 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:2000; ICC/IF-1:200-1:500; FC-1:10-1:100

Target information (Q15465)

Gene Name	SHH
About Gene Symbol	SHH (Sonic Hedgehog Signaling Molecule) is a metabolic enzyme involved in Hedgehog signaling. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Sonic Hedgehog in HepG2, Hela, 3T3 lysates using Sonic Hedgehog antibody.