

SAFB1 Rabbit mAb


ABZ2092

Product Overview


Name	SAFB1 Rabbit mAb
Accession(Primary)	Q15424
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 103 kDa; Observed MW: 150 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500; IF-1:50-1:100

Target information (Q15424)


Gene Name	SAFB
About Gene Symbol	



Immunocytochemistry analysis of SAFB1 (green)
in 293 using SAFB1 antibody, and DAPI (blue).



Western blot analysis of SAFB in K562, 3T3, HeLa
lysates using SAFB antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer tissue using SAFB antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.