

SF3B3 Rabbit mAb


ABZ2094

Product Overview


Name	SF3B3 Rabbit mAb
Accession(Primary)	Q15393
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 136 kDa; Observed MW: 136 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:500-1:2000; ICC/IF-1:200-1:500; FC-1:100-1:200; IP-1:20-1:50

Target information (Q15393)


Gene Name	SF3B3
About Gene Symbol	



Western blot analysis of SF3B3 in Raw264.7
lysates using SF3B3 antibody.



Western blot analysis of SF3B3 in SH-SY5Y, C6
lysates using SF3B3 antibody.



Immunocytochemistry analysis of SF3B3(green) in Hela using SF3B3 antibody, and DAPI(blue)