


ANGPT1 Rabbit mAb

ABZ2095

Product Overview

Name	ANGPT1 Rabbit mAb
Accession(Primary)	Q15389
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 58 kDa; Observed MW: 75 kDa
Recommended Dilution	WB-1:2000-1:5000

Target information (Q15389)

Gene Name ANGPT1**About Gene Symbol**

Western blot analysis of Angiopoietin 1 in Jurkat, C6, HeLa lysates using Angiopoietin 1 antibody.

