

TOMM20 Rabbit mAb


ABZ2096

Product Overview


Name	TOMM20 Rabbit mAb
Accession(Primary)	Q15388
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat, Hamster
Calculated MW	Calculated MW: 16 kDa; Observed MW: 16 kDa
Recommended Dilution	WB-1:1000-1:10000; IHC-1:100-1:500; IF: 1:50; FC-1:20-1:50

Target information (Q15388)


Gene Name	TOMM20
About Gene Symbol	



Western blot analysis of TOMM20 in HeLa, CHO-K1 lysates using TOMM20 antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using TOMM20 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of TOMM2 (green)
in HeLa using TOMM20 antibody, and DAPI(blue)