

PCBP2/hnRNP E2 Rabbit mAb


ABZ2097

Product Overview

Name	PCBP2/hnRNP E2 Rabbit mAb
Accession(Primary)	Q15366
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 39 kDa; Observed MW: 35-45 kDa
Recommended Dilution	WB-1:500-1:1000; ICC/IF-1:200-1:500; FC-1:100-1:200; IP-1:50-1:100

Target information (Q15366)

Gene Name	PCBP2
About Gene Symbol	PCBP2 (Poly(RC) Binding Protein 2) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of PCBP2/hnRNP E2 in C6, 3T3 lysates using PCBP2 antibody.