

IRF4 Rabbit mAb


ABZ2100

Product Overview

Name	IRF4 Rabbit mAb
Accession(Primary)	Q15306
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 52 kDa; Observed MW: 52 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; IP-1:10-1:100

Target information (Q15306)

Gene Name	IRF4
About Gene Symbol	



Western blot analysis of MUM1 in Ramos lysates
using Interferon Regulatory Factor 4 antibody.