

ErbB 4 Rabbit mAb

ABZ2101


Product Overview

Name	ErbB 4 Rabbit mAb
Accession(Primary)	Q15303
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 147 kDa; Observed MW: 180 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:200-1:1000

Target information (Q15303)

Gene Name ERBB4

About Gene Symbol



Western blot analysis of ErbB 4 in HeLa, A549, HL-60, U2OS, C6 lysates using ErbB 4 antibody.

