

nmt55/p54nrb Rabbit mAb


ABZ2102

Product Overview

Name	nmt55/p54nrb Rabbit mAb
Accession(Primary)	Q15233
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 54 kDa; Observed MW: 54 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500; ICC/IF-1:500-1:2000

Target information (Q15233)

Gene Name	NONO
About Gene Symbol	



Western blot analysis of NMT55/p54nrb in 293T, MOLT4 lysates using NONO antibody.