

p23 Rabbit mAb


ABZ2103

Product Overview


Name	p23 Rabbit mAb
Accession(Primary)	Q15185
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 19 kDa; Observed MW: 23 kDa
Recommended Dilution	WB-1:2000-1:5000; ICC/IF-1:100-1:200

Target information (Q15185)


Gene Name	PTGES3
About Gene Symbol	PTGES3 (Prostaglandin E Synthase 3) is a biosynthetic enzyme involved in receptor-mediated signaling. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Immunocytochemistry analysis of p23 (green) in Hela using p23 antibody, and DAPI (blue).



Western blot analysis of p23/PTGES3 in mouse brain, rat brain, mouse lung lysates using p23/PTGES3 antibody.



Western blot analysis of p23/PTGES3 in HeLa, A549, HL-60 lysates using p23/PTGES3 antibody